FavorPrep™ 96-Well Genomic DNA Kit

Cat: FADWE 96001 (1 回分) / FADWE 96002 (2 回分) / FADWE 96004 (4 回分)

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





v 221116

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	FADWE 96001 (1 prep)	FADWE 96002 (2 preps)	FADWE 96004 (4 preps)
FATG1 Buffer	40 ml	80 ml	80 ml × 2
FATG2 Buffer	40 ml	80 ml	80 ml × 2
W1 Buffer (concentrated)*	44 ml	88 ml	88 ml × 2
Wash Buffer (concentrated)*	17.5 ml	35 ml	35 ml × 2
Elution Buffer	30 ml	60 ml	60 ml × 2
Proteinase K**	24 mg	48 mg	48 mg × 2
Filter Plate (96-Well DNA Binding Plate)	1 plate	2 plates	4 plates
Collection Plate (96-Well 2ml Plate)	3 plates	6 plates	12 plates
Elution Plate (96-Well PCR Plate)	1 plate	2 plates	4 plates
Adhesive Film	2 pcs	4 pcs	8 pcs
*添加する 96-100%エタノール量			
W1 Buffer (concentrated)	16 ml	32 ml	32 ml × 2
Wash Buffer (concentrated)	70 ml	140 ml	140 ml × 2
*添加する ddH₂O 量			
Proteinase K	2.4 ml	4.8 ml	4.8 ml × 2

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<p>• STEP 1. Sample preparation and lysis</p> <p>For whole blood, buffy coat, serum, plasma, body fluids</p> <ul style="list-style-type: none"> • Add Proteinase K to a Collection Plate (first collection plate) • Add sample • Add FATG2 Buffer and mix completely • Seal with Adhesive Film. • Incubate the plate with shaking at 60°C for 20 min 					
<p>For aminam tissues and cultured cells</p> <ul style="list-style-type: none"> • Collect samples in a Collection Plate (first collection plate) • Add FATG1 Buffer and Proteinase K • Seal with Adhesive Film. Incubate the plate with shaking at 60°C for 1~2 hrs • Add FATG2 Buffer and mix completely • Seal with Adhesive Film. Further incubate the plate at 70°C for 20 min 					
<p>• STEP 2. Adjust binding condition:</p> <ul style="list-style-type: none"> • Add ethanol • Mix by pipetting 					
<p>• STEP 3. Bind DNA to Filter Plate:</p> <table border="1"> <thead> <tr> <th>centrifuge protocol</th> <th>vacuum protocol</th> </tr> </thead> <tbody> <tr> <td> <ul style="list-style-type: none"> • Combine the plates. • Transfer the sample mixture to Filter plate. • Centrifuge at 4,500 ~ 6,000 x g for 5 min.  </td> <td> <ul style="list-style-type: none"> • Fix Plates to manifold. • Transfer the sample mixture to Filter plate. • Apply 10 inches Hg vacuum until the wells have emptied.  </td> </tr> </tbody> </table>		centrifuge protocol	vacuum protocol	<ul style="list-style-type: none"> • Combine the plates. • Transfer the sample mixture to Filter plate. • Centrifuge at 4,500 ~ 6,000 x g for 5 min. 	<ul style="list-style-type: none"> • Fix Plates to manifold. • Transfer the sample mixture to Filter plate. • Apply 10 inches Hg vacuum until the wells have emptied. 
centrifuge protocol	vacuum protocol				
<ul style="list-style-type: none"> • Combine the plates. • Transfer the sample mixture to Filter plate. • Centrifuge at 4,500 ~ 6,000 x g for 5 min. 	<ul style="list-style-type: none"> • Fix Plates to manifold. • Transfer the sample mixture to Filter plate. • Apply 10 inches Hg vacuum until the wells have emptied. 				
<p>• STEP 4. Wash the Filter Plate Twice (W1 Buffer and Wash Buffer)</p> <table border="1"> <tbody> <tr> <td> <ul style="list-style-type: none"> • Add W1 Buffer. Centrifuge at 5,600 ~ 6,000 x g for 2 min • Add Wash Buffer. Centrifuge at 5,600 ~ 6,000 x g for 15 min  </td> <td> <ul style="list-style-type: none"> • Add W1 Buffer. Apply vacuum at 10 inches Hg. • Add Wash Buffer. Apply vacuum at 10 inches Hg for 10 min.  </td> </tr> </tbody> </table>		<ul style="list-style-type: none"> • Add W1 Buffer. Centrifuge at 5,600 ~ 6,000 x g for 2 min • Add Wash Buffer. Centrifuge at 5,600 ~ 6,000 x g for 15 min 	<ul style="list-style-type: none"> • Add W1 Buffer. Apply vacuum at 10 inches Hg. • Add Wash Buffer. Apply vacuum at 10 inches Hg for 10 min. 		
<ul style="list-style-type: none"> • Add W1 Buffer. Centrifuge at 5,600 ~ 6,000 x g for 2 min • Add Wash Buffer. Centrifuge at 5,600 ~ 6,000 x g for 15 min 	<ul style="list-style-type: none"> • Add W1 Buffer. Apply vacuum at 10 inches Hg. • Add Wash Buffer. Apply vacuum at 10 inches Hg for 10 min. 				
<p>• STEP 5. Dry the membranes of Filter Plate:</p> <table border="1"> <tbody> <tr> <td> <ul style="list-style-type: none"> • Stand the Filter plate on a clean paper towel at room temperature for 10 min. </td> <td> <ul style="list-style-type: none"> • Tap the Filter Plate tips on paper towel • Return the Filter Plate and the Collection Plate to the manifold. • Apply maximum vacuum for an additional 10 min. </td> </tr> </tbody> </table>		<ul style="list-style-type: none"> • Stand the Filter plate on a clean paper towel at room temperature for 10 min. 	<ul style="list-style-type: none"> • Tap the Filter Plate tips on paper towel • Return the Filter Plate and the Collection Plate to the manifold. • Apply maximum vacuum for an additional 10 min. 		
<ul style="list-style-type: none"> • Stand the Filter plate on a clean paper towel at room temperature for 10 min. 	<ul style="list-style-type: none"> • Tap the Filter Plate tips on paper towel • Return the Filter Plate and the Collection Plate to the manifold. • Apply maximum vacuum for an additional 10 min. 				
<p>• STEP 6. DNA Elution:</p> <table border="1"> <tbody> <tr> <td> <ul style="list-style-type: none"> • Add Elution Buffer or ddH₂O to the Filter Plate. Stand for 3 min. • Centrifuge to elute DNA.  </td> <td> <ul style="list-style-type: none"> • Add Elution Buffer or ddH₂O to the Filter Plate. Stand for 3 min. • Close the manifold valve. Turn on the vacuum source to build up a vacuum to 15 inches Hg. • Open the manifold valve to apply vacuum to elute DNA. <p>Alternative: If the consistent volume of elutes are recommend use centrifuge protocol to proceed this elution step. (Page 3, STEP 6)</p>  </td> </tr> </tbody> </table>		<ul style="list-style-type: none"> • Add Elution Buffer or ddH₂O to the Filter Plate. Stand for 3 min. • Centrifuge to elute DNA. 	<ul style="list-style-type: none"> • Add Elution Buffer or ddH₂O to the Filter Plate. Stand for 3 min. • Close the manifold valve. Turn on the vacuum source to build up a vacuum to 15 inches Hg. • Open the manifold valve to apply vacuum to elute DNA. <p>Alternative: If the consistent volume of elutes are recommend use centrifuge protocol to proceed this elution step. (Page 3, STEP 6)</p> 		
<ul style="list-style-type: none"> • Add Elution Buffer or ddH₂O to the Filter Plate. Stand for 3 min. • Centrifuge to elute DNA. 	<ul style="list-style-type: none"> • Add Elution Buffer or ddH₂O to the Filter Plate. Stand for 3 min. • Close the manifold valve. Turn on the vacuum source to build up a vacuum to 15 inches Hg. • Open the manifold valve to apply vacuum to elute DNA. <p>Alternative: If the consistent volume of elutes are recommend use centrifuge protocol to proceed this elution step. (Page 3, STEP 6)</p> 				

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