FavorPrep™ PCR Clean Up 96-Well Kit

Cat. No. : FAPC107A (1 回分) / FAPC107B (2 回分) / FAPC107C (4 回分)

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STEP 1. Sample preparation

Transfer samples in a Collection Plate (first Collection Plate)

Add FAPC Buffer

Mix by pipetting

Vacuum protocol


or

Centrifuge protocol

STEP 2. Bind DNA to Filter Plate:


Vacuum processing

- Transfer the sample mixture to Filter plate.
- Apply -12 inches Hg vacuum until the wells have emptied.

Centrifuge processing

- Transfer the sample mixture to Filter plate.
- Centrifuge at 4,500~6,000 xg for 2 mins.

STEP 3. Wash the Filter Plate with Wash Buffer

- Add Wash Buffer. Apply vacuum at -12 inches for 2 mins.

- Add Wash Buffer. Centrifuge at 5,600~6,000 x g for 10 mins.

Collection Plate (third Collection Plate)


STEP 4. Dry the membranes of the Filter Plate:

- Tap the Filter Plate tips on paper towel
- Return the Filter Plate and the second Collection Plate back to the manifold.
- Apply vacuum at -12 inches Hg for an additional 10 mins.


- Stand the Filter plate on a clean paper towel at room temperature for 5 min.

STEP 5. DNA Elution:

- Add Elution Buffer to the Filter Plate. Stand for 3 mins.
 - Close the manifold valve. Turn on the vacuum source to build up a vacuum to -12 inches Hg.
 - Open the manifold valve to apply vacuum to elute DNA.
- Alternative:** If the consistent volume of elutes are recommended, use centrifuge processing for this elution step. (Page 3, STEP 6)

- Add Elution Buffer to the Filter Plate. Stand for 3 mins.
- Centrifuge to elute DNA.

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