

Sample amount and typical yield

Sample	Recommended amount of sample used	Typical yield (µg)
Animal cells (up to 1×10 ⁷)	HeLa, 1×10 ⁶ cells	10
High yield Tissue (Mouse)(up to 20 mg)	Liver, 10 mg Spleen, 10 mg	35 45
Low yield Tissue (Mouse)(up to 50 mg)	Embryo, 10 mg Heart, 10 mg Brain, 10 mg Kidney, 10 mg Lung, 10 mg Intestine, 10 mg	10 7.5 7.5 20 10 15

Safety Information:

1. FARB Buffer and Wash Buffer 1 provided in this system contain irritants. Wear gloves and lab coat when handling these buffers.
2. **CAUTION:** FARB Buffers and Wash Buffer1 contain guanidinium salts which can form highly reactive compounds when combined with bleach. **DO NOT add bleach or acidic solutions directly to the preparation waste.**

Kit Component: FARB Buffer	
Hazard contents Guanidinium thiocyanate CAS-No. 593-84-0 EC-No. 209-812-1	
	 Danger
Hazard statement(s) H302 + H312 + H332 H314 H412	Harmful if swallowed, in contact with skin or if inhaled. Causes severe skin burns and eye damage. Harmful to aquatic life with long lasting effects.
Precautionary statement(s) P260 P280 P301 + P312 + P330 P303 + P361 + P353 P304 + P340 + P310 P305 + P351 + P338	Do not breathe dust/ fume/ gas/ mist/ vapours/ spray. Wear protective gloves/ protective clothing/ eye protection/ face protection. IF SWALLOWED: Call a POISON CENTER/ doctor if you feel unwell. Rinse mouth. IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. IF INHALED: Remove person to fresh air and keep comfortable for breathing. Immediately call a POISON CENTER/ doctor. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Kit Component: Wash Buffer 1	
Hazard contents Guanidine hydrochloride CAS-No. 50-01-1 EC-No. 200-002-3	
	 Warning
Hazard statement(s) H302 + H332 H315 H319	Harmful if swallowed or if inhaled. Causes skin irritation. Causes serious eye irritation.
Precautionary statement(s) P261 P301 + P312 + P330 P305 + P351 + P338	Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray. IF SWALLOWED: Call a POISON CENTER/ doctor if you feel unwell. Rinse mouth. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Protocol: Vacuum processing

Please read Important Notes and Additional Materials Required before starting the following steps.

STEP 1. Sample preparation and lysis

- For animal cells:**
- Transfer up to 1×10⁷ cells to each well of a Collection Plate. (provided, 96-well 2 ml plate; first Collection Plate). Centrifuge the plate at 500 ×g, 4°C for 5 mins. Remove the supernatant.
 - Add 450 µl of FARB Buffer (β-Me added). Pipet up and down to resuspend the cells completely.
 - Incubate the sample mixture at room temperature for 5 mins.
- For animal tissues:**
- Transfer up to 50 mg tissue to each well of a Collection Plate. (provided, 96-well 2 ml plate; first Collection Plate).
 - Add 450 µl of FARB Buffer (β-Me added).
 - Disrupt the sample with a appropriate homogenizer.
 - Incubate the sample mixture at room temperature for 5 min.

STEP 2. Clarify lysate

- Seal with the Adhesive Film on the Collection Plate. Place the plate in a rotor bucket and centrifuge at 5,600~6,000 ×g for 10 mins.

STEP 3. Adjust binding condition

- Transfer 350 µl of the upper clarified lysate to each well of a clean Collection Plate (provided, second Collection Plate).
-- **Note: Avoid to pipet any debris and pellet when transferring the supernatant.**
- Add 350 µl of RNase-free 70% ethanol to each well and mix by pipetting.
-- **Note: make sure that ethanol mixed with lysate completely.**

STEP 4. RNA Binding

- Fix a clean Collection Plate (provided, third Collection Plate) on the rack of vacuum manifold and cover the manifold lid. Place a Filter Plate (provided, 96-Well RNA binding plate) on top of the third Collection Plate.
- Transfer the sample mixture to the Filter Plate and discard the second Collection Plate).
- Apply vacuum at -12 inches Hg until the wells have emptied.
- Release vacuum from the manifold.
- Discard the Collection Plate (third).
- Place the Filter Plate and a clean Collection Plate (provided, fourth Collection Plate) to the manifold.

(Optional STEP): Digest DNA by DNase I

- Follow the steps from A1~A4 to eliminate DNA. Otherwise, proceed STEP 5 directly.**
- A1.** Add 250 µl of Wash Buffer 1 (ethanol added) to each well of the Filter Plate. Apply vacuum at -12 inches Hg for 2 mins.
Release vacuum from the manifold.
 - A2.** Add 750 µl of RNase-free 70% ethanol to each well of the Filter Plate. Apply vacuum at -12 inches Hg for 2 mins.
Release vacuum from the manifold.
 - A3.** Add 60 µl of DNase I reaction mixture (0.25 U/µl, not provided) to each well's membrane of the Filter Plate. Stand the plate for 15 mins at room temperature. **Do not vacuum after incubation. Proceed step A4 directly.**
 - A4.** Add 250 µl of Wash Buffer 1 to each well of the Filter Plate. Apply vacuum at -12 inches Hg until the wells have emptied. Release vacuum from the manifold. Discard the flow-through. Return the Filter Plate and Collection Plate back to the manifold.
 - A5. After DNase I treatment, proceed STEP 6.**

STEP 5. Wash the Filter Plate with Wash Buffer 1

- Add 500 µl of Wash Buffer 1 (ethanol added) to each well of the Filter Plate.
- Apply vacuum at -12 inches Hg until the wells have emptied.
- Release vacuum from the manifold.
- Discard the flow-through. Return the Filter Plate and the Collection Plate back to the manifold.

STEP 6. Wash the Filter Plate with Wash Buffer 2

- Add 500 µl of Wash Buffer 2 (ethanol added) to each well of the Filter Plate.
- Apply vacuum at -12 inches Hg for 2 mins.

- Release vacuum from the manifold.
- Discard the flow-through. Return the Filter Plate and the Collection Plate back to the manifold.

STEP 7. Wash the Filter Plate “again” with Wash Buffer 2

- Repeat Step 6

STEP 8. Dry the membranes of Filter Plate

- Gently tap the tips of the Filter Plate on a clean paper towel to remove residual liquid.
- Return the Filter Plate to the Collection Plate fixed in the manifold.
- Apply vacuum at -12 inches Hg for an addition **10 mins.**
- Release vacuum from the manifold.
- Discard the Collection Plate.

STEP 9. RNA Elution

- Alternative: If the consistent volume of eluates are recommended, use “centrifuge processing step 9-1~9-4”, to proceed this elution.**
- Place an Elution Plate (provided, 96-Well PCR Plate) on the 96-Well PCR Rack (not provided) and fix plate onto manifold. Cover the manifold lid and place the Filter Plate on the Elution Plate (top: Filter Plate; middle: 96-Well PCR Plate; bottom: 96-Well PCR Rack).
 - Add 50~75 µl of RNase-Free Water to the membrane center of the Filter Plate. Stand for 3 mins.
-- **Note! The eluates averaged about 25 µl less than the adding volume of elution buffers. For example, adding 50 µl of RNase-free water will recover ~25 µl of eluate.**
-- **Note! Do not use RNase-Free Water less than the suggested volume (<50 µl). It will lower the RNA yield.**
-- **Note! For effective elution, make sure that RNase-free water is dispensed on the membrane center and is absorbed completely.**
 - Close the manifold valve. Turn on the vacuum source to build up a vacuum to -12 inches Hg.
 - Open the manifold volve to apply vacuum to elute RNA.
 - Release vacuum from the manifold.
 - Take out the Elution Plate (96-well PCR plate) and seal with an Adhesive Film (provided). Store the RNA at -70°C before use.

Protocol: Centrifuge processing

Please read Important Notes and Additional Materials Required before starting the following steps.

STEP 1. Sample preparation and lysis

- For animal cells:**
- Transfer up to 1×10⁷ cells to each well of a Collection Plate. (provided, 96-well 2 ml plate; first Collection Plate). Centrifuge the plate at 500 ×g, 4°C for 5 mins. Remove the supernatant.
 - Add 450 µl of FARB Buffer (β-Me added). Pipet up and down to resuspend the cells completely.
 - Incubate the sample mixture at room temperature for 5 mins.
- For animal tissues:**
- Transfer up to 50 mg tissue to each well of a Collection Plate. (provided, 96-well 2 ml plate; first Collection Plate).
 - Add 450 µl of FARB Buffer (β-Me added).
 - Disrupt the sample with a appropriate homogenizer.
 - Incubate the sample mixture at room temperature for 5 mins.

STEP 2. Clarify lysate

- Seal with the Adhesive Film on the Collection Plate. Place the plate in a rotor bucket and centrifuge at 5,600~6,000 ×g for 10 mins.

STEP 3. Adjust binding condition

- Transfer 350 µl of the upper clarified lysate to each well of a clean Collection Plate (provided, second Collection Plate).
-- **Note: Avoid to pipet any debris and pellet when transferring the supernatant.**
- Add 350 µl of RNase-free 70% ethanol to each well and mix by pipetting.
-- **Note: make sure that ethanol mixed with lysate completely.**

STEP 4. RNA Binding

- Place a Filter Plate (provided, 96-Well RNA binding plate) on a clean Collection Plate (provided, third Collection Plate).

- Transfer the sample mixture to each well of the Filter Plate and discard the second Collection Plate.
- Place the combined plates (Filter Plate + the third Collection Plate) in a rotor bucket and centrifuge at 5,600~6,000 ×g for 2 mins.
- Discard the Collection Plate.
- Place the Filter Plate on a clean Collection Plate (provided, fourth Collection Plate).

(Optional STEP): Digest DNA by DNase I
Follow the steps from B1~B4 to eliminate DNA. Otherwise, proceed STEP 5 directly.

- B1.** Add 250 µl of Wash Buffer 1 (ethanol added) to each well of the Filter Plate. Place the combined plates in a rotor bucket and centrifuge at 5,600~6,000 ×g for 5 mins. Discard the flow-through and return the Filter Plate back to the Collection Plate.
- B2.** Add 750 µl of RNase-free 70% ethanol to each well of the Filter Plate. Place the combined plates in a rotor bucket and centrifuge at 5,600~6,000 ×g for 5 mins. Discard the flow-through and return the Filter Plate back to the Collection Plate.
- B3.** Add 60 µl of DNase I solution (0.25 U/µl, not provided) to each well's membrane of the Filter Plate. Stand the plates for 15 mins at room temperature. **Do not centrifuge after incubation. Proceed step B4 directly.**
- B4.** Add 250 µl of Wash Buffer 1 to each well of the Filter Plate. Place the plates in a rotor bucket and centrifuge at 5,600~6,000 ×g for 2 mins. Discard the flow-through and return the Filter Plate to the Collection Plate.
- B5. After DNase I treatment, proceed STEP 6.**

STEP 5. Wash the Filter Plate with Wash Buffer 1

- Add 500 µl of Wash Buffer 1 (ethanol added) to each well of the Filter Plate.
- Place the combined plate in a rotor bucket and centrifuge at 5,600~6,000 ×g for 2 mins.
- Discard the flow-through and return the Filter Plate back to the Collection Plate.

STEP 6. Wash the Filter Plate with Wash Buffer 2

- Add 500 µl of Wash Buffer 2 (ethanol added) to each well of the Filter Plate.
- Place the combined plate in a rotor bucket and centrifuge at 5,600~6,000 ×g for 2 mins.
- Discard the flow-through and return the Filter Plate to the Collection Plate.

STEP 7. Wash the Filter Plate “again” with Wash Buffer 2

- Add 500 µl of Wash Buffer 2 (ethanol added) to each well of the Filter Plate.
- Place the combined plate in a rotor bucket and centrifuge at 5,600~6,000 ×g for **10 mins.**
- Discard the Collection Plate.

STEP 8. Dry the membranes of Filter Plate

- Place the Filter Plate on top of a clean paper towel (not provided) and stand at room temperature for **5 mins.**

STEP 9. RNA Elution

- 9-1. Place the combined Filter Plate and Elution Plate (provided, 96-Well PCR Plate) onto the 96-Well PCR Rack (not provided), forming a three-plate assembly in the following order: top – DNA Binding Plate; middle – Elution Plate; bottom – 96-Well PCR Rack.
- 9-2. Add 50~75 µl of RNase-Free Water to the membrane center of the Filter Plate. Stand for 3 mins.
-- **Note! The eluates averaged about 25 µl less than the adding volume of elution buffers. For example, adding 50 µl of RNase-Free Water will recover ~ 25 µl of eluate.**
-- **Note! Do not use RNase-free water less than the suggested volume (<50 µl). It will lower the RNA yield.**
-- **Note! For effective elution, make sure that RNase-free water is dispensed on the membrane center and is absorbed completely.**
- 9-3. Place the combined plates in a rotor bucket and centrifuge at 5,600~6,000 ×g for 5 mins to elute RNA.
- 9-4. Take out the Elution Plate (96-well PCR plate) and seal with an Adhesive Film (provided). Store the RNA at -70°C before use.