

2023年3月28日

お客様各位

株式会社チヨダサイエンス

フェーバージェン社 核酸抽出製品  
構成品変更のご案内

拝啓 時下益々ご清栄のこととお慶び申し上げます。  
平素は格別のご高配を賜り、厚く御礼申し上げます。

さて、下記の通りFAVORGEN (フェーバージェン)社「核酸抽出キット」の一部製品につきまして、構成品が変更となりますので、ご案内申し上げます。

また、これに伴う価格の変更はございません。新しい取扱説明書は弊社WEBサイトにて公開しております。

大変ご迷惑をおかけいたしますが、ご理解と変わらぬご愛顧を賜りますようお願い申し上げます。

敬具

記

1. 対象製品

FAVORGEN (フェーバージェン) 「核酸抽出キット」一部製品  
別紙一覧表を参照ください。

2. 変更内容

Proteinase K (lyophilized → liquid)

3. 変更時期

2023年4月3日(月)以降入荷分より順次

以上

[お問い合わせ先]

営業企画部

TEL: 03-5294-7701 FAX: 03-5294-7752

E-mail: technical@chiyoda-s.jp

URL: <http://www.chiyoda-s.jp/>

## 【別紙】 FAVORGEN 構成品変更製品 一覧

品名	型番
Blood Genomic DNA Extraction Mini Kit (50)	FABGK 001
Blood Genomic DNA Extraction Mini Kit (100)	FABGK 001-1
Blood Genomic DNA Extraction Mini Kit (300)	FABGK 001-2
Blood Genomic DNA Extraction Maxi Kit (10)	FABGK 003
Blood Genomic DNA Extraction Maxi Kit (24)	FABGK 003-1
Tissue Genomic DNA Extraction Mini Kit (50)	FATGK 001
Tissue Genomic DNA Extraction Mini Kit (100)	FATGK 001-1
Tissue Genomic DNA Extraction Mini Kit (300)	FATGK 001-2
Tissue Genomic DNA Extraction Micro Elute Kit (50)	FATGM 001
Tissue Genomic DNA Extraction Micro Elute Kit (100)	FATGM 001-1
Stool DNA Isolation Kit (50)	FASTI 001
Stool DNA Isolation Kit (100)	FASTI 001-1
96-well Genomic DNA Extraction Kit (1 plate)	FADWE 96001
96-well Genomic DNA Extraction Kit (2 plates)	FADWE 96002
96-well Genomic DNA Extraction Kit (4 plates)	FADWE 96004